Tallas	Referencia proveedor	Referencia fabricante	código GTIN/EAN
5,5	GLE.101.01	85-173-1	9555076108330
6	GLE.102.01	85-173-2	9555076108347
6,5	GLE.103.01	85-173-3	9555076108354
7	GLE.104.01	85-173-4	9555076108361
7,5	GLE.105.01	85-173-5	9555076108378
8	GLE.106.01	85-173-6	9555076108385
8,5	GLE.107.01	85-173-7	9555076108392
9	GLE.108.01	85-173-8	9555076108408

Características

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MAXITEX®

PREPOWDERED SURGICAL GLOVES
STERILE

This product is classed as Personal Protective Equipment (PPE) according to European PPE Regulation 2016/425 and has shown to comply with this Regulation through the Harmonised standards EN ISO 374-1:2016, EN ISO 374-5:2016, EN 420:2003+A1:2009

PERFORMANCE AND LIMITATION OF USE - This product has been tested in accordance with EN ISO 374-1:2016, EN ISO 374-5:2016, EN 374-4:2013 and EN 420:2003+A1:2009 and achieved the following performance levels:

EN ISO 374-1:2016/Type B

Code Letter	Chemical	Permeation Level
K	Sodium hydroxide 40%	6
P	Hydrogen peroxide 30%	6
T	Formaldehyde 37%	6

EN ISO 374-5:2016

Test acc. to EN 374-5:2016 Level 2 (AQL 0.65) Water leak tested by manufacturer
Protection against bacteria and fungi - Pass
Protection against viruses - Pass

Performance level	AQL	Inspection levels
Level 3	<0.65	G1
Level 2	<1.5	G1
Level 1	<4.0	S4

Notified Body responsible for certification and for module C2 compliance is SATRA,

2777

SATRA Technology Centre
Bracelown Business Park
Clonee, D15 YWZP, Ireland,
NB# 2777.

EN ISO 374-1:2016 Permeation level are based on breakthrough times as follows:
Performance level 1 2 3 4 5 6
Measured breakthrough times (mins) >10 >30 >60 >120 >240 >480

EN 421:2010 for protection against particulate radioactive contamination

CONCLUSION:

When assessed in accordance with EN 374-4:2013 the samples of gloves described as Latex surgical gloves powdered (Maxitex PP) achieved the following degradation results:

Chemical	Mean degradation / %
40% Sodium hydroxide (CAS: 1310-73-2)	-17.0
30% Hydrogen peroxide (CAS: 7722-84-1)	-1.4
37% Formaldehyde (CAS: 50-00-0)	-0.9

NOTE: The quoted mean degradation values are subject to often significant measurement uncertainties. Please see results tables below for more information.

Tests acc. to EN 374-4:2013. Protective gloves against chemicals and micro-organisms, Part 4: Determination of resistance to degradation by chemicals.

RESULTS:

Sample description:	Latex surgical gloves powdered Maxitex PP								
	40% Sodium hydroxide			30% Hydrogen peroxide			37% Formaldehyde		
Challenge chemical:	(23 ± 1)			(23 ± 1)			(23 ± 1)		
Test temperature / °C:	Glove 1	Glove 2	Glove 3	Glove 1	Glove 2	Glove 3	Glove 1	Glove 2	Glove 3
Degradation / %:	-14.8	-13.3	-22.9	-1.1	-3.5	-0.1	-8.8	-2.7	1.8
Mean degradation (DR) / %:	-17.0			-0.9			-1.4		
Standard deviation (σ _m) / %:	5.1			0.7			6.4		
UoM * / ± %:	21.6			24.6			19.1		
Appearance of samples after testing:	No change			Swollen and discoloured			Swollen and discoloured		

Absolute measurement uncertainty of the mean degradation value; it is therefore inferred that the true degradation value, with 95% confidence, lies within the range (DR ± UoM) %.

NOTE: Where the test specimens gave an increased puncture force after chemical exposure, the result is reported as a negative degradation.

EN374-4:2013 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

Declaration of conformity can be found at www.adventa-health.com/ppe/eudoc/doc_maxitex.pdf

Available in sizes 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0