Tallas	Referencia proveedor	código GTIN/EAN
XS (5-6)	GNE.000.10	9555076101649
S (6-7)	GNE.111.10	9555076101652
M (7-8)	GNE.222.10	9555076101669
L (8-9)	GNE.333.10	9555076101670
XL (9-10)	GNE.444.10	9555076101683

Características

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**NUGARD® NITRIL
STERILE**

www.adventa-health.com

Nitrile examination gloves and protective product equipment of complex design. Powder-free, non-sterile, ambidextrous. The gloves fit for special purpose in accordance to the MDD EN 455-2, as the product is an examination glove, as well as for use for hand protection when handling chemicals. These gloves will not provide protection in the wrist area against chemical splashes. These gloves contain no substances known to cause allergies. A list of ingredients is available on request. Usage and Care: Gloves should be stored in a Cool, dry environment out of direct sunlight. Avoid contact with petroleum products and grease. Gloves life depends upon storage and wearing conditions. The user should make a visual inspection before using. Discard if damaged. They should not be washed or chemically cleaned without specialist advice. Gloves should only be transported in their original packaging. Gloves meeting the requirements of EN ISO 374-1:2016 and EN ISO 374-5:2016. The gloves are designed to protect against chemical risks in accordance with EN ISO 374-1:2016, to protect against micro-organisms risks in accordance with EN ISO 374-1:2016 general requirements in accordance with EN420:2003+A1:2009. Certification in accordance with PPE Regulation (EU) 2016/425. Individuals with a proven or suspected hypersensitivity to natural rubber latex or Type IV Sensitivity to Chemicals should consult a physician before using any products. Performance is tested on the palm of the hand and performance levels are shown as explained below. Available sizes: XS=6, S=7, M=8, L=9, XL=10.

EN ISO 374-1:2016/Type B

Tests acc. to EN ISO 374-1:2016

Code Letter	Chemical	Permeation Level
K	40% Sodium hydroxide	6
O	25% Ammonium hydroxide	0
P	30% Hydrogen peroxide	3
T	37% Formaldehyde	5

Permeation Performance level	Measured breakthrough time (mins)
1	>10
2	>30
3	>60
4	>120
5	>240
6	>480

EN ISO 374-5:2016

Tested in accordance to EN 374-5:2016
Protection against bacteria and fungi - PASS
Protection against viruses - PASS
Performance level: LEVEL 2 (ISO 2859, AQL 1.0, G1)

Performance level	AQL	Inspection levels
Level 3	<0,65	G1
Level 2	<1,5	G1
Level 1	<4,0	S4

Compliant with:

MEDICAL	MDD 93/42/EEC, 2007/47/EC	EN455 - 1,2,3,4	Class I
PPE	(EU)2016/425	EN 420:2003+A1:2009 EN ISO 374-1:2016 EN ISO 374-5:2016	Category III

Compliant with ASTM F1671

EN 421:2010

Notified Body responsible for certification and for module C2 compliance is SATRA.

2777

SATRA Technology Centre
Bracelown Business Park
Clonee, D15 YN2P, Ireland.
NB# 2777.

Terang Nusa Sdn Bhd
1 Jalan 8
Pengkalan Chepa 2 Industrial Zone
16100 Kota Bharu, Malaysia

EC REP ulma International GmbH
Pfaffenweg 35
89231 Neu-Ulm, Germany
Made in Malaysia

CONCLUSION:

When assessed in accordance with EN 374-4:2013 the samples of gloves described as ENRBF50-Light blue nitrile gloves achieved the following degradation results:

Chemical	Mean degradation / %
40% Sodium hydroxide (CAS: 1310-73-2)	-17,0
30% Hydrogen peroxide (CAS: 7722-84-1)	-1,4
37% Formaldehyde (CAS: 50-00-0)	-0,9

NOTE: The quoted mean degradation values are subject to often significant measurement uncertainties. Please see results tables below for more information.

Tests acc. to EN 374-4:2013. Protective gloves against chemicals and micro-organisms. Part 4: Determination of resistance to degradation by chemicals.

RESULTS:

Sample description:	ENRBF50-Light blue											
	40% Sodium hydroxide			25% Ammonium hydroxide			30% Hydrogen peroxide			37% Formaldehyde		
Challenge chemical:	(23 ± 1)			(23 ± 1)			(23 ± 1)			(23 ± 1)		
Test temperature / °C:												
Degradation / %:	Glove 1	Glove 2	Glove 3	Glove 1	Glove 2	Glove 3	Glove 1	Glove 2	Glove 3	Glove 1	Glove 2	Glove 3
	-16.6	-2.9	6.9	40.5	40.0	44.7	21.8	42.7	34.0	27.0	37.2	26.6
Mean degradation (DR) / %:	-4.2			41.8			32.9			30.3		
Standard deviation (σ _{DR}) / %:	11.8			2.6			10.5			6.0		
UoM * / ± %:	22.1			27.3			18.8			18.4		
Appearance of samples after testing:	Swollen			Swollen			No change			Swollen		

Absolute measurement uncertainty of the mean degradation value; it is therefore inferred that the true degradation value, with 95% confidence, lies within the range (DR ± UoM) %.

NOTE: Where the test specimens gave an increased puncture force after chemical exposure, the result is reported as a negative degradation.

EN374-4:2013 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

Declaration of conformity can be found at:
www.adventa-health.com/ppe/eu/doc/doc_nugard_nitril_sterile.pdf